

# OTImix

## Automatic system for synthesis of Nitrox



OTImix is a last-generation automatic system for synthesis of Nitrox breathing mixtures [ $N_2 + O_2$ ] used in hyperbaric oxygen-therapy applications.

Differently from traditional solutions, OTImix employs two proportional solenoid valves managed by a microcontroller, so that  $O_2$  concentration and pressure in mixture remain aligned to the set points, constantly and independently from the outlet flow.

New solutions allows replacing the classical lunging sequence by a continue flow from the mixing circuit, coinciding with the current average consumption.

Advantages of this organization are small size of storing tank, very good accuracy, stability and repeatability in the synthesis process, direct use of different diluent gases as Air,  $N_2$  and He.

OTImix implements two different startup sequences, based on partial pressure calculation, in order to supply a fast correction of the remaining mixture in tank.

The **ECOLOGICAL** startup sequence saves all the contents of the tank, compensating the current consumption by  $O_2$  or Air /  $N_2$ , until the  $O_2$  concentration in the tank reaches the current set point.


The **STANDARD** startup sequence partially discharges the tank, so that the remaining mixture can reach the desired concentration and pressure through only one fast charge by  $O_2$  or Air /  $N_2$ .

In case the initial conditions of the remaining mixture allow reaching the desired conditions without discharging the tank, both sequences execute one or two fast charges of  $O_2$  and/or Air /  $N_2$  only, before starting the operating status.

User can set all the reference and system parameters by five pushbuttons on the control panel. The only manual operation is the adjusting of the inlet pressure regulator.

The *Knudsen* effect  $O_2$  analyzer, integrated in the equipment, verifies the composition of the mixture and supplies excellent performances regarding accuracy, long-term stability, operating temperature, independence from atmospheric pressure and flow, response time. It needs a periodic maintenance.

## TECHNICAL CHARACTERISTICS



*OTImix control panel*

- Min P inlet (Air / N<sub>2</sub> / He/ O<sub>2</sub>)
- Max P inlet (Air / N<sub>2</sub> / He/ O<sub>2</sub>)
- P outlet (°)
- Mix O<sub>2</sub> concentration [% vol] (°)
- Mix composition combined error
- Max outlet flow

P outlet + 1 bar  
15 bar  
from 5.5 to 9 bar  
from 25 to 90 %  
< ± 1% O<sub>2</sub> vol  
150 NI/m O<sub>2</sub> +  
150 NI/m N<sub>2</sub> /Aria

(°) – Programmable by user

- Board size
 


|     |         |
|-----|---------|
| L = | 500 mm  |
| P = | 350 mm  |
| H = | 1200 mm |
- Proof degree
 

|  |       |
|--|-------|
|  | IP 22 |
|--|-------|
- Tank capacity
 

|  |      |
|--|------|
|  | 50 l |
|--|------|
- Supply
 

|  |         |
|--|---------|
|  | 220 Vac |
|--|---------|
- Max power
 

|  |      |
|--|------|
|  | 50 W |
|--|------|



*OTImix connections*